



# VaporFit User Guide

**Przemysław Pastwa and Piotr Bruździak**

Department of Physical Chemistry, Gdańsk University of Technology,  
Narutowicza 11-12, 80-233 Gdańsk, Poland  
piotr.bruzdziak@pg.edu.pl

## Contents

<b>1</b>	<b>Introduction</b>	<b>2</b>
<b>2</b>	<b>License</b>	<b>2</b>
<b>3</b>	<b>Citations</b>	<b>2</b>
<b>4</b>	<b>Data sets</b>	<b>2</b>
<b>5</b>	<b>User Guide</b>	<b>3</b>
5.1	VaporFit window overview . . . . .	3
5.2	Load Spectral Files . . . . .	4
5.3	Load Atmospheric Spectra . . . . .	4
5.4	Set SG Parameters and Wavenumber Range . . . . .	4
5.5	Run the Correction Algorithm . . . . .	4
5.6	(Optional) Run “Optimum Search” Tool . . . . .	5
5.6.1	Smoothness Functions . . . . .	5
5.7	(Optional) Perform PCA Analysis . . . . .	7
5.8	(Optional) Check Correction Coefficients . . . . .	8
<b>6</b>	<b>Save Results</b>	<b>9</b>
6.1	Save Report . . . . .	9
6.2	Save Corrected Spectra . . . . .	9
<b>7</b>	<b>Exit the Program</b>	<b>9</b>
<b>8</b>	<b>Final Notes</b>	<b>10</b>

## 1 Introduction

VaporFit is an open-source software tool designed for automated atmospheric correction in FTIR spectroscopy. It employs a **least squares fitting approach** using multiple reference atmospheric/vapor spectra, dynamically adjusting subtraction coefficients to improve accuracy. You are welcome to send us any comments or suggestions.

## 2 License

This software is distributed under the **GNU GPL v.3.0** license + citation requirement (see LICENSE file). By using this software, you agree to its terms and conditions.

## 3 Citations

If you use VaporFit in your research, you must cite the following publications:

- Bruzdziak P., Vapor correction of FTIR spectra – A simple automatic least squares approach, *Spectrochimica Acta Part A: Molecular and Biomolecular Spectroscopy*, 2019, **223**, 1–4, doi:10.1016/j.saa.2019.117373
- Pastwa P., Bruzdziak P., VaporFit: An Open-Source Software for Accurate Atmospheric Correction of FTIR Spectra, *in preparation*

## 4 Data sets

Three sample data sets are available to download from the project site:

- **test\_data\_betaine\_set.zip** — a set of 29 spectra of betaine aqueous solutions with two atmospheric spectra ( $2\text{ cm}^{-1}$ , 128 scans, single background measurement),
- **test\_data\_D2O-H2O\_set.zip** — a set of 11 spectra of  $\text{D}_2\text{O}/\text{H}_2\text{O}$  mixtures with three atmospheric spectra ( $4\text{ cm}^{-1}$ , 128 scans, single background measurement).
- **test\_data\_urea\_evaporation\_set.zip** – a set of 20 spectra of a  $20\text{ }\mu\text{L}$  urea solution (initial concentration  $1\text{ M}$ ), freely evaporated on the ATR crystal for 20 minutes, with two atmospheric spectra (resolution:  $4\text{ cm}^{-1}$ , 64 scans, single background measurement).

All spectra were measured on a Bruker Invenio-R FTIR spectrometer equipped with a single reflection diamond ATR accessory. The spectrometer was purged with dry nitrogen (Claind Brezza NiGen LCMS 40.1) for 30 minutes before and during the measurement.

## 5 User Guide

### 5.1 VaporFit window overview




Figure 1: **Left:** tools and parameters. **Right:** original spectra before correction (upper left panel), atmospheric spectra (upper right panel), correction results (bottom panel; solid lines—corrected spectra, dashed lines—spectra before corrections).

Use the **left mouse button** to select and zoom a fragment of the spectrum (each panel is independent). Use the **right mouse button** to reset the view. Zooming in or out does not affect the wavenumber range for correction indicated in the section of the VaporFit window.

## 5.2 Load Spectral Files

Upload spectra that require correction.

**Supported formats:**

- CSV
- DPT

**Important:** Ensure that the spectral resolution, wavenumber range, number of points, and units are **consistent across all files** to avoid errors.

## 5.3 Load Atmospheric Spectra

Upload atmospheric spectra.

- The format and spectral parameters must match those of the spectra being corrected.

## 5.4 Set SG Parameters and Wavenumber Range

Adjust the **desired wavenumber range** and the **Savitzky-Golay (SG) parameters** in the Parameters Section of the main window.


- Default settings are predefined for standard atmospheric corrections for the amide I region of the spectra. Adjust as needed.
- If values exceed the range of loaded spectra, they will be limited to the closest minimum or maximum wavenumber.
- It is **not recommended** to perform correction over the entire spectrum range (typically 4000–400  $\text{cm}^{-1}$ ). Choosing a narrower region improves result quality.
- Modifying SG parameters may improve results, but incorrect settings can introduce artifacts. A good rule of thumb is to use the lowest possible SG parameter values for optimal correction.

## 5.5 Run the Correction Algorithm

Press “**Correct!**” to apply the atmospheric correction algorithm.

- The resulting graph will display:
  - **Original spectra** (black dashed lines)
  - **Corrected spectra** (color-coded)
- The correction process **automatically optimizes subtraction coefficients** based on least squares fitting.

## 5.6 (Optional) Run “Optimum Search” Tool



This tool generates a **graph showing the dependence of four spectral smoothness metrics on the SG window size (+/-10 points) for a given SG polynomial order**. Values represent the sum of these metrics for all corrected spectra.

- The minima of these functions (shown in the legend) suggest **potential optimal SG window sizes**.
- However, these values **should not be used blindly**. Verify with **PCA analysis** to ensure meaningful correction.
- The four available smoothness metrics may have different applications depending on the spectral characteristics. For amide I band spectra, the SD parameter works best.
- Depending on computational power, this step may take some time.

### 5.6.1 Smoothness Functions

#### 1. Standard Deviation of the Residue After Smoothing (SD):

$$\sigma_{\text{residue}} = \text{std}(\text{spectrum} - \text{smoothed spectrum})$$

where smoothed spectrum is the result of applying the Savitzky-Golay filter to the spectrum Spectrum (with default values of polynomial order 3 and window width 11).

## 2. Spectral Smoothness Index (SSI):

$$\text{SSI} = \frac{\sum \left(\frac{d}{dx}\text{spectrum}\right)^2}{\sum (\text{spectrum})^2}$$


where  $\sum$  denotes the sum of all data points of the spectrum or its function, and  $\frac{d}{dx}\text{spectrum}$  is the first derivative of the spectrum, i.e., the differences between consecutive values of the spectrum.

## 3. Second Derivative Variance (SDV):

$$\text{Variance of second derivative} = \text{Var} \left( \frac{d^2}{dx^2}\text{spectrum} \right)$$

where  $\frac{d}{dx}\text{spectrum}$  is the second derivative of the spectrum, i.e., the differences between consecutive values of the  $\frac{d}{dx}\text{spectrum}$ .

## 5.7 (Optional) Perform PCA Analysis



After correction, press **“PCA”** to analyze the **principal components (PCs)** of the spectral series **before** (left) and **after** (right) correction. You can set the number of extracted PC spectra in the main window. Use the **left mouse button** to select and zoom a fragment of the spectrum. Use the **right mouse button** to reset the view. Zooming in or out in one spectrum graph will affect the other.


**PCA window includes:**

- **Explained Variance (EV) plots** (top)—A sharp drop suggests the **minimum number of significant PCs**.
- **Principal component shapes** (bottom)—Before correction, atmospheric components should be recognizable in PCs. After proper correction, they should no longer be present.

**Important Notes:**

- PCA **will not run** if the correction is incomplete or if fewer than **two spectra** are available.
- Misinterpretation of PCA can lead to incorrect assumptions about correction quality.
- PCA in this software is for visual assessment only. The actual number and shape of principal components should be determined using specialized software or appropriate libraries.

## 5.8 (Optional) Check Correction Coefficients



Use the **“Parameters Graph”** tool to inspect the correction coefficients applied to each atmospheric reference spectrum.



## 6 Save Results

### 6.1 Save Report

Press **“Save report”** to export **multiple CSV files**, which include:

- General information about the correction (date, number of spectra, file names, and paths)
- Spectra **before and after correction** (as a matrix; for the selected wavenumber range)
- **Atmospheric reference spectra** within the same range (as a matrix)
- **Correction parameters** for each spectrum
- **PCA results**, if performed

You will be asked for the report filename. The given name will serve as a prefix to all these files.

### 6.2 Save Corrected Spectra

Press **“Export Spectra”** to save all corrected spectra, each in a separate CSV file, in the same folder as the original spectra. All filenames will start with the “corr\_” prefix.

## 7 Exit the Program

Press **“Close”** to exit VaporFit.

## 8 Final Notes

- VaporFit is designed for atmospheric corrections in **FTIR spectroscopy** but can process any CSV spectral files.
- If you experience issues with correction quality, double-check **input spectra consistency** and consider adjusting **wavenumber ranges** or **reference spectra selection**.
- If correction is unsatisfactory, ensure you have enough atmospheric spectra covering varied conditions; set the SG polynomial order to 3 and the SG window size to 5, or try correction with a single atmospheric spectrum.
- The software is actively developed, and new features may be added in future releases.
- We tested and compiled the software on MacOS v.15.3.1 (Mac M1 M4) and Windows 11 (Intel Core i5-8365).
- Try running the Python code in your Python environment if you are experiencing issues with compiled versions.